

Triple Mix Formulations we have compounded at Partners In Care

Papaverine 12mg, Phentolamine 0.4 mg, Prostaglandin 4mcg/ml	<input type="checkbox"/> 5ml	<input type="checkbox"/> 10 ml
Papaverine 16.7 mg, Phentolamine 0.55mg, Prostaglandin 5.67mcg/ml	<input type="checkbox"/> 5ml	<input type="checkbox"/> 10 ml
Papaverine 17.65mg, Phentolamine 0.59mg, Prostaglandin 6mcg/ml	<input type="checkbox"/> 5ml	<input type="checkbox"/> 10 ml
Papaverine 17.79mg, Phentolamine 0.62 mg, Prostaglandin 6mcg/ml	<input type="checkbox"/> 5ml	<input type="checkbox"/> 10 ml
Papaverine 18mg, Phentolamine 1mg, Prostaglandin 6mcg/ml	<input type="checkbox"/> 5ml	<input type="checkbox"/> 10 ml
Papaverine 18mg, Phentolamine 1mg, Prostaglandin 8mcg/ml	<input type="checkbox"/> 5ml	<input type="checkbox"/> 10 ml
Papaverine 22.5mg, Phentolamine 0.83 mg, Prostaglandin 8.33mcg/ml	<input type="checkbox"/> 5ml	<input type="checkbox"/> 10 ml
Papaverine 24mg, Phentolamine 0.9 mg, Prostaglandin 10mcg/ml	<input type="checkbox"/> 5ml	<input type="checkbox"/> 10 ml
Papaverine 27mg, Phentolamine 1mg, Prostaglandin 20mcg/ml	<input type="checkbox"/> 5ml	<input type="checkbox"/> 10 ml
Papaverine 30mg, Phentolamine 0.2mg, Prostaglandin 10 mcg/ml	<input type="checkbox"/> 5ml	<input type="checkbox"/> 10 ml
Papaverine 30mg, Phentolamine 0.5 mg, Prostaglandin 10 mcg/ml	<input type="checkbox"/> 5ml	<input type="checkbox"/> 10 ml
Papaverine 30mg, Phentolamine 0.5mg, Prostaglandin 20mcg/ml	<input type="checkbox"/> 5ml	<input type="checkbox"/> 10 ml
Papaverine 30mg, Phentolamine 1 mg, Prostaglandin 2mcg/ml	<input type="checkbox"/> 5ml	<input type="checkbox"/> 10 ml
Papaverine 30mg, Phentolamine 1mg, Prostaglandin 10mcg/ml	<input type="checkbox"/> 5ml	<input type="checkbox"/> 10 ml
Papaverine 30mg, Phentolamine 2mg, Prostaglandin 20mcg/ml	<input type="checkbox"/> 5ml	<input type="checkbox"/> 10 ml
Papaverine 30mg, Phentolamine 1.5mg, Prostaglandin 50mcg/ml	<input type="checkbox"/> 5ml	<input type="checkbox"/> 10 ml
Papaverine mg, Phentolamine mg, Prostaglandin mcg/ml	<input type="checkbox"/> 5ml	<input type="checkbox"/> 10 ml